

50LD

LN₂ DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	50LD
LN ₂ Capacity (L)	50
Static Hold Time (Days)	122
Evap. Rate (Liters/Day)	0.41
DIMENSIONS	
Necktube (in/mm)	2.5 / 63.5
Overall Dia. (in/mm)	18.8 / 477.5
Height (in/mm)	32.4 / 823
WEIGHT	
Empty (lb/kg)	38.7 / 17.6
Full (lb/kg)	127.7 / 57.9

FEATURES

- 50 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency