

5LD

LN₂ DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	5LD
LN ₂ Capacity (L)	5
Static Hold Time (Days)	6
Evap. Rate (Liters/Day)	0.77
DIMENSIONS	
Necktube (in/mm)	5.6 / 142.2
Overall Dia. (in/mm)	7.6 / 193
Height (in/mm)	17.5 / 444.5
WEIGHT	
Empty (lb/kg)	6.9 / 3.1
Full (lb/kg)	15.8 / 7.2

FEATURES

- 5 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency