

# LS750

## LABORATORY SYSTEMS

For liquid or vapor phase vial storage in convenient box-type racks.



SPECIFICATIONS	LS750
LN <sub>2</sub> Capacity (L)	35
Static Hold Time (Days)	130
Evap. Rate (Liters/ day)	0.27
DIMENSIONS	
Necktube (in/mm)	4.7 / 119.4
Overall Dia. (in/mm)	18.8 / 477.5
Height (in/mm)	26.8 / 680.7
WEIGHT	
Empty (lb/kg)	39 / 17.7
Full (lb/kg)	101.3 / 45.9

## FEATURES

- Convenient Storage Systems – rack index location ring and internal location fixture, computer compatible box storage for simple inventory management
  - Number of Racks: 6
  - Number of Shelves / Rack: 5
  - 2 ml Vial Capacity: 750
- In vapor-phase storage, temperature at the top of racks is -190°C.
- High strength aluminum body and magneformed necktube design
- Superior vacuum performance with super insulation provides maximum holding times
- Lockable lids