

# 10XT

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	10XT
LN <sub>2</sub> Capacity (L)	10
Static Hold Time (Days)	100
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	11.4 / 289.6
Height (in/mm)	23.8 / 604.5
WEIGHT	
Empty (lb/kg)	16.5 / 7.5
Full Weight (lb/kg)	34.3 / 15.6
CAPACITY	
Number of Canisters	6
Canister Size (in/mm)	1.5 x 11 / 38 x 279
1.2 & 2.0 ml Vials (5/cane)	150
1.2 & 2.0 ml Vials (6/cane)	180
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	N/A
.5 cc Straws – Bulk (2 Levels)	1500



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)

[www.ICBiomedical.com](http://www.ICBiomedical.com)

Made In USA

© 2022, IC Biomedical, LLC