

# 20XT

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	20XT
LN <sub>2</sub> Capacity (L)	20.7
Static Hold Time (Days)	230
Evap. Rate (Liters/Day)	0.09
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	15.6 / 396.2
Height (in/mm)	25.8 / 655.3
WEIGHT	
Empty (lb/kg)	26 / 11.8
Full Weight (lb/kg)	62.8 / 28.5
CAPACITY	
Number of Canisters	6
Canister Size (in/mm)	1.5 x 11 / 38 x 279
1.2 & 2.0 ml Vials (5/cane)	150
1.2 & 2.0 ml Vials (6/cane)	180
.5 cc Straws on Canes	540
.5 cc Straws – Bulk (1 Level)	750
.5 cc Straws – Bulk (2 Levels)	1500