

3XTL

LN₂ EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	3XTL
LN ₂ Capacity (L)	3
Static Hold Time (Days)	27
Evap. Rate (Liters/Day)	0.11
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	7.6 / 193
Height (in/mm)	17.2 / 436.8
WEIGHT	
Empty (lb/kg)	7.2 / 3.2
Full Weight (lb/kg)	12.5 / 5.7
CAPACITY	
Number of Canisters	6
Canister Size (in/mm)	1.5 x 11 / 38 x 279
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A