

# 8XTL

## LN<sub>2</sub> EXTENDED TIME SMALL FREEZER

For long-term storage at cryogenic temperatures.



### FEATURES

- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Versatile storage system
- Superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	8XTL
LN <sub>2</sub> Capacity (L)	8
Static Hold Time (Days)	80
Evap. Rate (Liters/Day)	0.1
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	15.6 / 396.2
Height (in/mm)	19 / 482.6
WEIGHT	
Empty (lb/kg)	19.6 / 8.9
Full Weight (lb/kg)	33.8 / 15.3
CAPACITY	
Number of Canisters	6
Canister Size (in/mm)	1.5 x 11 / 38 x 279
.5 cc Straws on Canes	N/A
.5 cc Straws – Bulk (1 Level)	750
Number of 1.2 & 2.0 ml Vials	N/A