

CLASSIC-25

LN₂ DEWAR

For safe storage and dispensing of small accounts of liquid nitrogen.



SPECIFICATIONS	CLASSIC-25
LN ₂ Capacity (L)	25
Static Hold Time (Days)	119
Evap. Rate (Liters/Day)	0.21
DIMENSIONS	
Necktube (in/mm)	2.0 / 50.8
Overall Dia. (in/mm)	15.5 / 393.7
Height (in/mm)	22.9 / 581.6
WEIGHT	
Empty (lb/kg)	19 / 8.6
Full (lb/kg)	63.5 / 28.8

FEATURES

- 25 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com

Made In USA

© 2022, IC Biomedical, LLC