

CXR 100

VAPOR (DRY) SHIPPER

For the safe transport of biological samples in a vapor phase.



FEATURES

- 3.7 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

SPECIFICATIONS	CXR100
LN ₂ Capacity (L)	3.7
Static Hold Time (Days)	16
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in/mm)	3.58 / 90.9
Overall Dia. (in/mm)	9.2 / 233.7
Height (in/mm)	19.4 / 492.7
WEIGHT	
Empty (lb/kg)	11.7 / 5.3
Charged Weight (lb/kg)	18.2 / 8.3
CAPACITY	
No. of 1.2 ml & 2.0 ml vials (5/cane)	85
No. of 1.2 ml & 2.0 ml vials (6/cane)	102
Number of .25 cc Straws (28/cane)	784
No. of .25 cc Straws, Bulk (2 levels) in Goblets	1640
Number of 1/2 straws (10/cane)	280
No. of .5 cc Straws, Bulk (2 levels) in Goblets	800
Number of Canisters	1