

CXR500

VAPOR (DRY) SHIPPER

For the safe transport of biological samples in a vapor phase.



FEATURES

- 6.4 Liter absorbed LN2 capacity
- Advanced Concept Adsorbent enables faster charging
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid

| SPECIFICATIONS | CXR500 |
|--|--------------|
| LN ₂ Capacity (L) | 6.4 |
| Static Hold Time (Days) | 11 |
| Evap. Rate (Liters/Day) | 0.64 |
| DIMENSIONS | |
| Necktube (in/mm) | 8.5 / 216 |
| Overall Dia. (in/mm) | 15.5 / 393.7 |
| Height (in/mm) | 26.9 / 683.3 |
| WEIGHT | |
| Empty (lb/kg) | 30 / 13.6 |
| Charged Weight (lb/kg) | 41.5 / 18.8 |
| CAPACITY | |
| No. of 1.2 ml & 2.0 ml vials (5/cane) | 500 |
| No. of 1.2 ml & 2.0 ml vials (6/cane) | 500 |
| Number of .25 cc Straws (28/cane) | 1736 |
| No. of .25 cc Straws, Bulk (2 levels) in Goblets | 7160 |
| Number of 1/2 straws (10/cane) | 1240 |
| No. of .5 cc Straws, Bulk (2 levels) in Goblets | 3200 |
| Number of Canisters | N/A |