



Series

LN2 FREEZERS



Specifically designed for the every day requirements of laboratory and biorepository operations

- Universal canister compatibility
- Fits in standard shipping boxes and containers
- Rugged ribbed high-strength aluminum body, magneformed necktube design and more durable paint
- Superior vacuum performance with insulation technology that delivers maximum holding times
- Lockable lid with enlarged opening for seals or securing ties

| SPECIFICATIONS | TW20 | TW24 | TW34 | TW47-6 | TW47-10 |
|------------------------------|--------------|--------------|------------|--------------|--------------|
| LN ₂ Capacity (L) | 21 | 24 | 34 | 48 | 48 |
| Static Hold Time (Days) | 135 | 218 | 123 | 76 | 76 |
| Evap. Rate (Liters/Day) | 0.095 | 0.110 | 0.180 | 0.390 | 0.390 |
| DIMENSIONS | | | | | |
| Necktube (in/mm) | 2.2 / 55.9 | 2.38 (60.4) | 3.5 / 88.9 | 5 / 127 | 5 / 127 |
| Overall Dia. (in/mm) | 14.5 / 368.3 | 14.5 / 368.3 | 18 / 457.2 | 20 / 508 | 20 / 508 |
| Height (in/mm) | 26 / 660.4 | 27.3 / 693.4 | 27 / 685.8 | 26.5 / 673.1 | 26.5 / 673.1 |
| WEIGHT | | | | | |
| Empty (lb/kg) | 26 / 11.8 | 26.2 (11.9) | 35 / 15.9 | 31.5 / 14.3 | 31.5 / 14.3 |
| Full Weight (lb/kg) | 63 / 28.6 | 69.0 (31.3) | 94 / 42.6 | 120 / 54.4 | 120 / 54.4 |
| CAPACITY | | | | | |
| Number of Canisters | 6 | 6 | 6 | 6 | 10 |
| Canister Diameter (in/mm) | 1.65 / 41.9 | 1.89 / 48 | 2.8 / 71.2 | 4.0 / 101.6 | 2.8 / 71.1 |
| .5 cc Straws on Canes | 660 | | 2100 | 3900 | 3200 |

COMING SOON!

Model 47 dual-configuration for quick charge with leak back absorbent protection –avoid irritating hazardous shipment supplemental charges!



MADE IN USA



855.750.8191 | sales@ICBiomedical.com

www.ICBiomedical.com

Made In USA

© 2022, IC Biomedical, LLC