

Product Specifications

High-Capacity Laboratory Cryostorage Freezers

K-Series

DIMENSIONS		3K	10K	24K	38K
Overall Height	(in)	29.7	44.0	44.0	49.0
	(mm)	754	111.8	111.8	1245
Overall Dia.	(in)	15.4	23.1 x 30.5	34.0 x 38.5	42.0
	(mm)	391	587 x 775	864 x 965	1067
Neck Dia.	(in)	14.0	21.0	31.0	39.0
	(mm)	356	533	787	991
Internal Height	(in)	19.2	29.0	29.0	29.0
	(mm)	488	737	737	737
Internal Dia.	(in)	14.0	21.0	31.0	39.0
	(mm)	356	533	787	991



1-855-750-8191



www.icbiomedical.com



1031 Cass White Road, NW
Cartersville, Georgia 30121

Product Specifications

High-Capacity Cryostorage Freezers

LABS-Series

DIMENSIONS		20K	38K	40K	80K	80K-AI	AI-91K	94K
Neck Dia.	(in)	13	18	18	24.5	24.5	24.5	24.5
	(mm)	330	457	457	622	622	622	622
Overall Height	(in)	54	57	54	54	54	45	63
	(mm)	1371	1447	1371	1371	1371	1135	1600
Overall Dia.	(in)	34	42	44.6	59.5	59.5	60	59.5
	(mm)	864	1067	1143	1511	1511	1524	1511
Internal Height	(in)	30	30	30	30	30	13	35
	(mm)	762	762	762	762	762	332	897
Internal Dia.	(in)	29.5	37.5	40.5	55	55	56.5	55
	(mm)	750	965	1029	1397	1397	1443	1397
Step Height	(in)	11	11.5	11 / 16.3	11 / 17	-	-	11 / 22
	(mm)	275	287.5	275 / 407.5	275 / 425	-	-	275 / 550



1-855-750-8191



www.icbiomedical.com



1031 Cass White Road, NW
Cartersville, Georgia 30121

Product Specifications

High-Capacity Cryostorage Freezers

Precision & Revolution

DIMENSIONS		20K	24K	40K	45K	90K	100K
Neck Dia.	(in)	11.8	11.8	16.8	16.8	24.5	24.5
	(mm)	300	300	426	426	622	622
Overall Height	(in)	56.5	63	56	62	58	64
	(mm)	1437	1602	1419	1574	1469	1625
Overall Dia.	(in)	32	32	42	42	60	60
	(mm)	814	814	1068	1068	1524	1524
Useable Height	(in)	31	37	31	37	32	37.5
	(mm)	785	936	788	947	800	958
Internal Dia.	(in)	30	30	40	40	58	58
	(mm)	760	760	1012	1012	1470	1470



1-855-750-8191



www.icbiomedical.com



1031 Cass White Road, NW
Cartersville, Georgia 30121