

Product Specifications

Vapor (Dry) Shippers

CX/CXR & D/DX Series

| SPECIFICATIONS | | CX100 | CXR100 | CXR500 | 4D | 4DX | 10DX |
|---------------------|------|-------|--------|--------|-------|-------|-------|
| Evap. Rate (L/Day) | | 0.2 | 0.23 | 0.64 | 0.17 | 0.19 | 0.64 |
| Static Holding Days | | 18 | 16 | 11 | 27 | 24 | 11 |
| LN2 Capacity (L) | | 4.1 | 3.7 | 6.4 | 4.1 | 4.1 | 6.4 |
| DIMENSIONS | | | | | | | |
| Overall Height | (in) | 18.5 | 19.4 | 29.6 | 18.5 | 18.5 | 29.6 |
| | (mm) | 472.4 | 492.7 | 683.3 | 470 | 470 | 683.3 |
| Overall Dia. | (in) | 8.75 | 9.2 | 15.5 | 8.75 | 8.75 | 15.5 |
| | (mm) | 222.3 | 233.7 | 393.7 | 222.2 | 222.2 | 393.7 |
| Neck Dia. | (in) | 2.75 | 3.58 | 8.5 | 2.0 | 2.78 | 8.5 |
| | (mm) | 69.9 | 90.9 | 216 | 50.8 | 70.6 | 216 |



1-855-750-8191



www.icbiomedical.com



1031 Cass White Road, NW
Cartersville, Georgia 30121

Product Specifications

Vapor (Dry) Shippers

TW Dual Series

| SPECIFICATIONS | | Dual 20 | Dual 22 | Dual 34 | Dual 47-6 | Dual 47-10 |
|---------------------|------|---------|---------|---------|-----------|------------|
| Evap. Rate (L/Day) | | 0.095 | 0.180 | 0.180 | 0.390 | 0.390 |
| Static Holding Days | | 21 | 20 | 20 | 21 | 21 |
| LN2 Capacity (L) | | 21 | 21 | 33 | 48 | 48 |
| DIMENSIONS | | | | | | |
| Overall Height | (in) | 26 | 26 | 27 | 26.5 | 26.5 |
| | (mm) | 660.4 | 660.4 | 685.8 | 673.1 | 673.1 |
| Overall Dia. | (in) | 14.5 | 14.5 | 18 | 20 | 20 |
| | (mm) | 368.3 | 368.3 | 457.2 | 508 | 508 |
| Neck Dia. | (in) | 2.2 | 3.5 | 3.5 | 5.0 | 5.0 |
| | (mm) | 55.9 | 88.9 | 88.9 | 127 | 127 |



1-855-750-8191



www.icbiomedical.com



1031 Cass White Road, NW
Cartersville, Georgia 30121