

35LD

LN₂ DEWAR

For safe storage and dispensing of small amounts of liquid nitrogen.



SPECIFICATIONS	35LD
LN ₂ Capacity (L)	35
Static Hold Time (Days)	152
Evap. Rate (Liters/Day)	0.23
DIMENSIONS	
Necktube (in/mm)	2.5 / 63.5
Overall Dia. (in/mm)	18.8 / 477.5
Height (in/mm)	26.3 / 668
WEIGHT	
Empty (lb/kg)	35.1 / 15.9
Full (lb/kg)	97.4 / 44.2

FEATURES

- 35 Liter LN₂ capacity
- Rugged ribbed high-strength aluminum body, magneformed necktube design, and more durable paint
- Snap-on cap and necktube assures a tight closure
- Advanced insulation materials assure high thermal efficiency