

# LS Series

## VIAL STORAGE



For liquid- or vapor storage in convenient box type racks.

- Convenient Storage Systems – rack index location ring and internal location fixture, computer compatible box storage for simple inventory management
- In vapor-phase storage, temperature at the top of racks is  $-190^{\circ}\text{C}$ .
- High strength aluminum body and magneformed necktube design
- Superior vacuum performance with super insulation provides maximum holding times
- Lockable lids

SPECIFICATIONS	LS750	LS3000	LS4800	LS6000
LN <sub>2</sub> Capacity (L)	35	81	130	165
Static Hold Time (Days)	130	106	153	194
Evap. Rate (Liters/Day)	0.27	0.76	0.85	0.84
Number of Racks	6	6	6	6
Number of Shelves/Rack	5	5	8	10
2 ml Vial Capacity	750	3000	4800	6000
DIMENSIONS				
Necktube (in/mm)	4.7 / 119.4	8.5 / 215.9	8.5 / 215.9	8.5 / 215.9
Overall Dia. (in/mm)	18.8 / 477.5	26.9 / 683.3	26.9 / 683.3	26.9 / 683.3
Height (in/mm)	6.8 / 680.7	28.8 / 731.5	35.1 / 891.5	39 / 990.6
WEIGHT				
Empty (lb/kg)	39 / 17.7	70 / 31.8	90 / 40.8	121 / 54.9
Full Weight (lb/kg)	101.3 / 45.9	214.2 / 97.2	312.4 / 141.7	410 / 186



### Security Options

- Accessory low-level alarm is available on the LS 6000
- CS100 controller with remote connection capability is available on LS4800 and LS6000



MADE IN USA



855.750.8191 | [sales@ICBiomedical.com](mailto:sales@ICBiomedical.com)  
[www.ICBiomedical.com](http://www.ICBiomedical.com)

Made In USA

© 2023, IC Biomedical, LLC