

10DX

CRYO EXCHANGE VAPOR (DRY) SHIPPER

For the safe transport of biological samples in vapor phase.



SPECIFICATIONS	10DX
LN ₂ Capacity (L)	10
Static Hold Time (Days)	14
Evap. Rate (Liters/Day)	0.7
DIMENSIONS	
Necktube (in/mm)	8.5 / 215.9
Overall Dia. (in/mm)	15 / 381
Height (in/mm)	23.25 / 590.5
CAPACITY	
No. of 1.2 ml & 2.0 ml vials (6/cane)	966
Blood Bags	10

Note: Evaporation rate and static holding time are nominal. Actual rate may be affected by the nature of the contents, atmospheric conditions, container history and manufacturing tolerances.

FEATURES

- Hydrophobic Adsorbent material is easier to clean if the vessel gets contaminated
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid



855.750.8191 | sales@ICBiomedical.com
www.ICBiomedical.com
Made In USA