

# VAPOR SHIPPER UPRIGHT SHIPPING CASE



For the safe transport of vapor shippers

The new IC Biomedical Upright Shipping Case is made of high-impact-resistant material to protect its contents from physical damage. The Case's domed top prevents stacking and its circular wide base encourages upright orientation during handling and transit.

With shipping times getting longer and package handling getting rougher, IC Biomedical has designed this case to provide a safer, more secure enclosure for LN2 vapor shippers.

## Model 4LSC

For IC Biomedical vapor shipper models:

- 4D
- 4DX
- CX100
- CXR100

Height: 24-inches

Base Diameter: 15-inches

Empty Weight: 12 lbs

The 4LSC is compatible with most other similar-sized shippers available on the market.

## Model 10LSC

Designed for the 10DX vapor shipper.

Height: 28.3-inches

Base Diameter: 24-inches

Empty Weight: 19.98 lbs



Built-in slot allows the case to be locked down for added sample security.



800.886.2796 | Fax 317. 297.7988  
sales@ICBiomedical.com  
www.ICBiomedical.com