

D Series

VAPOR (DRY) SHIPPERS

For the safe transport of biological samples in a vapor phase.

- Hydrophobic Adsorbent material that's easier to clean if the vessel gets contaminated
- Complies with IATA regulations for open cryogenic receptacles
- Rugged construction and superior vacuum performance with super insulation for maximum holding times
- Lockable lid



SPECIFICATIONS	4D	4DX	10DX
LN ₂ Capacity (L)	3.7	3.7	10
Static Hold Time (Days)	21	19	14
Evap. Rate (Liters/Day)	0.17	0.19	0.7
DIMENSIONS			
Necktube (in/mm)	2.0/50.8	2.78/70.6	8.5/215.9
Overall Dia. (in/mm)	8.75/222.2	8.75/222.2	15/381
Height (in/mm)	18.5/470	18.5/470	23.25/590.5
WEIGHT*			
Empty (lb/kg)	11.7/5.31	10/4.53	24 / 10.9
Full Weight (lb/kg)	19/8.62	16.5/7.5	39.5 / 17.9
CAPACITY			
No. of 1.2 ml & 2.0 ml vials (5/cane)	40	85	N/A
No. of 1.2 ml & 2.0 ml vials (6/cane)	48	102	966
Number of .25 cc Straws (28/cane)	448	784	N/A
No. of .25 cc Straws, Bulk (2 levels) in Goblets	952	1640	N/A
Number of 1/2 straws (10/cane)	160	280	N/A
No. of .5 cc Straws, Bulk (2 levels) in Goblets	460	800	N/A
Number of Canisters/Blood Bags	1 Canister	1 Canister	10 Bags

* Weights are for guidance only. Actual weights may differ due to manufacturing and filling variances
 Note: Evaporation rate and static holding time are nominal. Actual rate may be affected by the nature of the contents, atmospheric conditions, container history and manufacturing tolerances.



MADE IN USA



855.750.8191 | sales@ICBiomedical.com | www.ICBiomedical.com

© 2023, IC Biomedical, LLC