



05 JANUARY 2023

TECHNICAL BULLETIN TB-2022-09-01 Revision B

SUBJECT: PRODUCTS SUPPLIED WITH ELECTRONIC CONTROLLERS - EUROPE

In an effort to improve customer delivery times and to streamline our product offering, with immediate effect IC Biomedical will now supply LABS-Series, K-Series and LS-Series with a Pacer manufactured control system.

For customers who still wish to use Mowden controllers and accessories which are purchased separately, we offer LABS and K-Series vessels that have the plumbing and solenoid valve ready for the installation of a Mowden M505CE controller.

LS series vessels will only be available with a CS100 controller going forward.

If you are interested in a -PL (plumbed) version of a vessel, item numbers are listed below and the distributor price lists will be updated to reflect the new item numbers.

10K-PL

24K-PL

38K-PL

LABS20K-PL

LABS38K-PL

LABS40K-PL

LABS80K-PL

LABS94K-PL

Please be aware that IC Biomedical can accept no liability for controllers purchased and installed by a third party, only for the vessel as supplied.

TECHNICAL SUPPORT

IC BIOMEDICAL, LLC