

August 22, 2023

TECHNICAL PRODUCT BULLETIN TB-2023-08-01 Rev A

Subject: 34HC Product Obsolescence

IC Biomedical would like to confirm the product obsolescence of the 34HC. The 34HC has now been replaced by the TW 34 vessel which has more durable plastics and an increased hold time.

The TW 34 has been specifically designed for the every day requirements of laboratory and biorepository operations.

Features include:

- *Universal canister compatibility*
- *Fits in standard shipping boxes and containers*
- *Rugged ribbed high-strength aluminum body, magnaformed necktube design and more durable paint.*
- *Superior vacuum performance with insulation technology that delivers maximum holding times.*
- *Lockable lid with enlarged opening seals for securing ties.*

SPECIFICATIONS	TW34
LN ₂ Capacity (L)	34
Static Hold Time (Days)	123
Evap. Rate (Liters/Day)	0.180
DIMENSIONS	
Necktube (in/mm)	3.5/88.9
Overall Dia. (in/mm)	18/457.2
Height (in/mm)	27/685.8
WEIGHT	
Empty (lb/kg)	35/15.9
Full Weight (lb/kg)	94/42.6
CAPACITY	
Number of Canisters	6
Canister Diameter (in/mm)	2.8/71.2

 1-855-750-8191

 www.icbiomedical.com

 1031 Cass White Road, NW
Cartersville, Georgia 30121