

# Safer Filling of Liquid Nitrogen Canisters for Cryotherapy Procedures

## Liquid Withdrawal Device (LWD)

The IC Biomedical Liquid Withdrawal Device (LWD) provides convenient, safer pressurized transfer of liquid nitrogen from a dewar to cryo gun canisters used in cryosurgical procedures.



Roller Base

## LWD FEATURES & BENEFITS

- Easily fill your handheld liquid nitrogen delivery system
- Permanently affixes to the top of the liquid nitrogen storage dewar and features a faucet-like on/off control to make the filling process simple
- Available for 25LD, 35LD and 50LD ICB Liquid Dewars
- Includes phase separator designed to minimize splashing and vaporization during filling
- Primary and secondary relief valves prevent hazardous pressure buildup.
- Easy filling for your LN2 supplier directly through the liquid outlet
- Safety ring and cable provides additional safety

Part #	Description
25LD	Liquid Dewar (25 liters)
35LD	Liquid Dewar (35 liters)
50LD	Liquid Dewar (50 liters)
D050-8C03	Liquid Withdrawal Device
R018-8C00	Roller Base 25LD
R033-8C00	Roller Base 35LD, 50LD

