



Product Specifications

Vapor (Dry) Shippers

CX/CXR & D/DX Series

SPECIFICATIONS		CX100	CXR100	CXR500	4D	4DX	10DX
Evap. Rate (L/Day)		0.2	0.23	0.64	0.17	0.19	0.64
Static Holding Days		18	16	11	27	24	11
LN2 Capacity (L)		3.7	3.7	6.4	4.1	4.1	6.4
DIMENSIONS							
Overall Height	(in)	18.5	19.4	29.6	18.5	18.5	29.6
	(mm)	472.4	492.7	683.3	470	470	683.3
Overall Dia.	(in)	8.75	9.2	15.5	8.75	8.75	15.5
	(mm)	222.3	233.7	393.7	222.2	222.2	393.7
Neck Dia.	(in)	2.75	3.58	8.5	2.0	2.78	8.5
	(mm)	69.9	90.9	216	50.8	70.6	216



Product Specifications Vapor (Dry) Shippers TW Dual Series

SPECIFICATIONS		Dual 20	Dual 22	Dual 34	Dual 47-6	Dual 47-10
Evap. Rate (L/Day)		0.095	0.180	0.180	0.390	0.390
Static Holding Days		21	20	20	21	21
LN2 Capacity (L)		21	21	33	48	48
DIMENSIONS						
Overall Height	(in)	26	26	27	26.5	26.5
	(mm)	660.4	660.4	685.8	673.1	673.1
Overall Dia.	(in)	14.5	14.5	18	20	20
	(mm)	368.3	368.3	457.2	508	508
Neck Dia.	(in)	2.2	3.5	3.5	5.0	5.0
	(mm)	55.9	88.9	88.9	127	127

Ref: PSD-002 Rev A



1-855-750-8191



www.icbiomedical.com



1031 Cass White Road, NW
Cartersville, Georgia 30121