



Initial Release January, 2024
Revised January, 2025

TECHNICAL PRODUCT BULLETIN TB-2024-01-02 Rev B

Subject: 35HC Product Obsolescence

IC Biomedical would like to confirm the product obsolescence of the 35HC. The 35HC has now been replaced by the ICB 35-10 vessel which has more durable plastics and an increased hold time.

The ICB 35-10 has been specifically designed for the every day requirements of laboratory and biorepository operations.

Features include:

- *Universal canister compatibility*
- *Fits in standard shipping boxes and containers*
- *Rugged ribbed high-strength aluminum body, magnaformed necktube design and more durable paint.*
- *Superior vacuum performance with insulation technology that delivers maximum holding times.*
- *Lockable lid with enlarged opening seals for securing ties.*

SPECIFICATIONS	ICB35-10
LN ₂ Capacity (L)	35
Static Hold Time (Days)	130
Evap. Rate (Liters/Day)	.027
DIMENSIONS	
Necktube (in/mm)	4.7 / 119.4
Overall Dia. (in/mm)	18 / 457.2
Height (in/mm)	27 / 685.8
WEIGHT	
Empty (lb/kg)	35 / 15.9
Full Weight (lb/kg)	96 / 43.5
CAPACITY	
Number of Canisters	10
Canister Diameter (in/mm)	2.64 x 67

 1-855-750-8191

 www.icbiomedical.com

 1031 Cass White Road, NW
Cartersville, Georgia 30121