



April, 2024

TECHNICAL PRODUCT BULLETIN TB-2024-04-01 Rev A

Subject: 38K Lid Design Change

In an effort to continuously improve the quality of the products IC Biomedical manufacture, we would like to confirm that there have been changes made to the 38K product. IC Biomedical have changed the lid material from plastic to metal, there is also the addition of sensor lines which will be used for future product development.

Specifications of the 38K have not changed, there are only visual changes.



 1-855-750-8191

 www.icbiomedical.com

 1031 Cass White Road, NW
Cartersville, Georgia 30121